





Revistin TM HEALBETTER, RECOVER QUICKER

What Is The ReviStim?



- It is a pulsed direct current stimulator
- It is in the E-stim categories of High Volt Pulsed Current (Galvanic) and Micro-Current devices

How Does It Work?

- It penetrates deep into the tissue, affecting cell and vessel
- It works below sensory

What Does The ReviStim Do?

- Increases Micro-circulation and new capillary formation
- Increases Oxygen levels in tissue
- Increases white blood cell migration
- Increases cellular production of ATP
- Kills bacteria
- Increases organization of cellular repair
- Increases fibroblasts migration to damaged site
- Increases strength of collagen deposits
- Reduces Edema and Inflammation
- Reduces Pain

How Easy is it to Use the ReviStim?



What the ReviStim Can Treat?

- Diabetic Neuropathy
- Muscle Pain
- Injuries: Whiplash, Sports, Work Comp
- Sprains/Strains
- Joint Pain
- Fractures: Slow healing, non-union
- Chronic Wounds/ Pressure Ulcers



Nocturnal Feature and Benefit of the ReviStim

- Because it works below sensory, it can be used while sleeping
- P2 mode allows for eight therapy sessions while sleeping

Why At Night?

- This allows for more therapy sessions
- This provides more efficient therapy because the body is more regenerative at night

Heal Better, Recover Quicker

